



Technical sheet
REVISION
 Ⓐ ARROW I



ARROW II option/revision tray allows insert and glenosphere replacements in case of ARROW I revision.


ARROW I revision
Keeping the humeral stem

REVERSE



ARROW II option/revision complement tray (*ref. A 271 907*)


ANATOMICAL



ARROW II option/revision complement tray + ARROW I universal anatomical tray (*ref. A 272 069*)


ARROW I revision
Keeping the glenoid part

REVERSE



ARROW II option/revision complement tray (*ref. A 271 907*)


ANATOMICAL



ARROW II option/revision complement tray
 + ARROW I universal anatomical tray (*ref. A 272 069*)



**ARROW I
IMPLANTS**



ARROW I UNIVERSAL ANATOMICAL TRAY
(ref. A 272 069)

ARROW I implants nomenclature
anatomical revision included
(ref. I 272 104)

ARROW II OPTION/REVISION TRAY
(ref. A 271 907)

ARROW II implants nomenclature
option/revision included
(ref. I 271 911)



ARROW I HUMERAL HEAD

Reference	Designation
260 532	ARROW centered humeral head D48 H18
260 535	ARROW centered humeral head D50 H19
260 550	ARROW centered humeral head D54 H19



ARROW I HUMERAL INSERT

Reference	Designation
257 057	ARROW humeral insert D36 H00
257 058	ARROW humeral insert D36 H05
257 059	ARROW humeral insert D36 H10
257 060	ARROW humeral insert D39 H00
257 061	ARROW humeral insert D39 H05
257 062	ARROW humeral insert D39 H10



ARROW I GLENOSPHERE

Reference	Designation
260 566	ARROW glenosphere D36
260 567	ARROW glenosphere D39



CENTERED HUMERAL HEAD

Reference	Diameter	Height (mm)
260 537	40	15
260 538	40	17
260 539	44	16
260 540	44	18
260 541	46	16
260 542	46	18
260 543	46	21
260 544	48	16
260 545	48	18
260 546	48	21
260 547	50	17
260 548	50	19
260 549	50	21
260 550	54	19
260 551	54	21



OFF-CENTERED HUMERAL HEAD

Reference	Diameter	Height (mm)
260 526	44	16
260 527	44	18
260 528	46	16
260 529	46	18
260 530	46	21
260 531	48	16
260 532	48	18
260 533	48	21
260 534	50	17
260 535	50	19
260 536	50	21



GLENOID INSERT

Reference	Size
260 556	44
260 557	46
260 558	48

**GROUPE
FH
ORTHO™**
AN OLYMPUS GROUP COMPANY



**DISTRIBUTEURS
DISTRIBUTORS**

FR, FH ORTHO SAS

3 rue de la Forêt - Zone Industrielle
BP 50009
68990 Heimsbrunn CEDEX - FRANCE
Tél. +33 (0)3 89 81 90 92
Fax : +33 (0)3 89 81 80 11
info@fhortho.com
www.fhortho.com

USA, FH ORTHO INC.

OrthoEx
16000 North 80th Street
Suite E
Scottsdale, Arizona 85260 - USA
Phone: +1 (412) 965-0950
customerservice@fhortho-us.com
www.fhortho.com

UK, FH ORTHO LTD

Unit 1b, Century Park, Valley way
Swansea Enterprise Park
Swansea, SA6 8RP - UK
Phone: +44 (0) 1792 464792
Fax: +44 (0) 844 412 7674
customer-servicesUK@fhortho.com
www.fhortho.com

PL, FH ORTHO POLSKA

Ul. Garbary 95/A6,
61-757 Poznan - POLSKA
Phone: +48 61 863 81 27
Fax: +48 61 863 81 28
biuro@implants24.pl
www.fhortho.com

FR, FH INDUSTRIE

6 rue Nobel, Z.I. de Kernevez
29000 QUIMPER - FRANCE
Tél. +33 (0)2 98 55 68 95
Fax : +33 (0)2 98 53 42 13
contact-fhi@fhortho.com
www.fhortho.com



**FABRICANT
MANUFACTURER**